

CAD/CAM DIGITAL SOLUTIONS

Digital Solutions




CAD/CAM Abutments

For all implant type diameters (3.2, 3.7, 4.2, 5.0mm)




Tray transfers are designed for impressions using the open-tray technique when the retentions are sharp, and for the closed-tray technique when they are round.

Material Titanium Alloy (Ti6Al4V ELI) & PEEK

Conical Connection

| | | Catalog Number |
|---|----------------------|----------------|
|  | Scan Body | Bio-SBC20 |
|  | Tabse Conical | Bio-TBC10 |
|  | Tabse no hex conical | Bio-TBN15 |

Internal Hex 2.42

| | | Catalog Number |
|---|--------------|----------------|
|  | Scan Body | Bio-SBH21 |
|  | Tabse hex | Bio-TBH11 |
|  | Tabse no hex | Bio-TBN16 |





info@biotec-implant.com
T. +49.6196.6521096
biotec-implant.com

Franzgraben 6,
34125 Kassel,
Germany

www.biotec-implant.com