

ONE PIECE ARP PROSTHETIC



One Piece | **ARP** | Prosthetics

Biotec offers a wide range of easy-to-use and esthetic restorations and prosthetics for almost any clinical need. The diameters of these prosthetics are interchangeable based on preference and need.



One Piece | ARP | Prosthetics

Transfers ARP

For all implant types diameters (3.2, 3.7, 4.2, 5.0mm)

Tray transfers are designed for impressions using the open-tray technique when the retentions are sharp, and for the closed-tray technique when they are round.



Material Body - Stainless Steel; Screws - Titanium

Transfer

Catalog Number



Bio-OPTRN

Drivers ARP

Catalog Number

Short



Bio-OPDRV10

Long



Bio-OPDRV15

Analog ARP

Implant Analogs (IAs) can be used for all implant diameters (3.2, 3.7, 4.2, 5.0mm)

we recommend using lab analog with identical dimensions to ensure the best rendering of the clinical situation.



Analog

Catalog Number



Bio-OPANG



info@biotec-implant.com
T. +49.6196.6521096
biotec-implant.com

Franzgraben 6,
34125 Kassel,
Germany

www.biotec-implant.com